


FRESE COMPATIBILI

COMPATIBLE DRILLS



CILINDRICHE CYLINDRICAL

Diametro Diameter	2	2,8	3,2	3,4	4,2	Altezza Height
CODICE CODE	FR20007	FR28007	FR32007	FR34007	FR42007	7
	FR20085	FR28085	FR32085	FR34085	FR42085	8,5
	FR20010	FR28010	FR32010	FR34010	FR42010	10
	FR20115	FR28115	FR32115	FR34115	FR42115	11,5
	FR20013	FR28013	FR32013	FR34013	FR42013	13
	FR20015	FR28015	FR32015	FR34015	FR42015	15

FISSE FIXED

Diametro Diameter	2	2,8	3,2	3,4	4,2	Altezza Height
CODICE CODE	FR200	FR280	FR320	FR340	FR420	-

- \varnothing 2 FRESA PILOTA
PILOT DRILL
- \varnothing 2,8 FRESA COMPATIBILE PER IMPIANTI \varnothing 3,25 Biomet 3i™
DRILL COMPATIBLE Biomet 3i™ IMPLANTS \varnothing 3,25
- \varnothing 3,2 FRESA COMPATIBILE PER IMPIANTI \varnothing 3,75 Biomet 3i™
DRILL COMPATIBLE Biomet 3i™ IMPLANTS \varnothing 3,75
- \varnothing 3,4 FRESA COMPATIBILE PER IMPIANTI \varnothing 4 Biomet 3i™
DRILL COMPATIBLE Biomet 3i™ IMPLANTS \varnothing 4
- \varnothing 4,2 FRESA COMPATIBILE PER IMPIANTI \varnothing 5 Biomet 3i™
DRILL COMPATIBLE Biomet 3i™ IMPLANTS \varnothing 5



IMPLANTOLOGIA SHOP

Viale della mercanzia, 113/115
40050 Funo, Argelato

e-mail: info@implantologiashop.it
tel.: +39 349 4095105