

BOX PER STRUMENTARIO CHIRURGICO

SURGICAL INSTRUMENT BOX

Box dal design ergonomico che si adatta perfettamente agli strumenti inseriti. Leggero, compatto e facile da trasportare.

Realizzato in plastica resistente agli urti, certificato per la sterilizzazione in autoclave.

Il tray porta-strumenti è stato realizzato per contenere accuratamente driver, strumentario chirurgico e frese, in sezioni separate.

Supporti in silicone individuali trattengono ogni strumento per evitarne il movimento durante il trasporto, anche se il box viene capovolto.

Le marcature laser sul tray sono anti-graffio, anti sbiadimento e consentono una facile pulitura.

MATERIALI

Box: resina polifenilsulfone Radel®

Porta-strumenti: silicone autoclavabile per uso medico

Vaschetta: acciaio inox

This ergonomically designed box is perfect for the tools it contains. It is light, compact and easy to move.

The box is made of crush resistant plastics and is autoclave-safe certified.

The tools tray was designed to contain drivers, surgical and milling tools, in separate sections. Individual silicone supports hold each tool, to prevent it from moving during transportation, even when the box is capsized.

The laser markings on the tray are scratch-proof and fade-proof and are easy to clean.

MATERIALS

Box: Radel® polyphenylsulfone

Support for instruments: autoclavable silicon for medical use

Metal bowl: stainless steel



CODICE BOX01
CODE



IMPLANTOLOGIA SHOP

Viale della mercanzia, 113/115
40050 Funo, Argelato

e-mail: info@implantologiashop.it
tel.: +39 349 4095105